Laboratoire

Laboratory

**Ref :
535 014**

Français – p 1

English – p 7

Version : 8007

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
** 0,15 € TTC/min. + prix appel*

1 Diagram

1. Case
2. Power Switch of Indicator
3. Temperature controller
4. Inner door knob
5. Toughened-glass inner door
6. Shelf
7. Magnetic air-tight loop of the outer door
8. Outer door

2 Scope of application

It is a necessary apparatus for storage of strains and culture of organisms in the fields of medical treatment, public health, pharmaceuticals, biology, agriculture and scientific research.

3 Technical features


Model	DNP-9022
Mains voltage	220 V, 50 Hz
Temperature range	Ambient temperature +5-65°C
Temperature fluctuation	± 0.5°C
Power consumption	<150 W
Dimensions of the working chamber (mm)	250 x 250 x 320

4 Structure


1. The incubator is electricity heated. The housing and the door are made of high quality steel plate and the outer surface of the unit is lacquer coated. The inner door was made of toughened glass. The culture samples in the chamber can be watched clearly without opening the inner glass door. The chamber was made of stainless steel. The stainless steel wire-mesh shelf is easy to move and the height is adjustable. A silicone gasket between the chamber and the inner glass door guarantees a satisfaction sealing. Three plate heaters were installed on the left, right and at the bottom at the outside of the chamber. A fan was installed on the top of the chamber running quietly, which also guarantee the uniformity of the inside atmosphere temperature. Superior glass fiber was filled between the chamber and outer casing for excellent insulation. There is an assistant heating system located at the bottom of the working chamber for special heating application. It can heat the chamber up to max. The power switch, power indicator and the microprocessor intelligent temperature controller were all installed on the top of the incubator. The microprocessor intelligent temperature controller is equipped with self-adjust PID method. It is of higher accuracy and rapid controlling, comparing to

normal PID controlling method. The set value of the temperature and the actual temperature are all digitally displayed. The controller can be easily to operate with soft-touch pad. Upper limit temperature limit as alarm function was included as standard function.

5 The function panel of the controller and the operation procedure

Operation procedure

Remark : Except the setting of temperature, all the other data of setting have been locked. If you wish to revise the new data.

6 Operation

1. Turn the power switch to '1' . The indicator of the power supplier will be light and there will be digital temperature display on the controller.

2. Temperature setting

After switched on the unit, if the set temperature is the same as the working temperature demanded, then it is not necessary to set it again. Otherwise, it is necessary to re-set the working temperature. First, press the function key 'SET' on the control panel to enter the temperature setting program, till the setting temperature display SV is flashing. Then, press the 'SHIFTING KEY' and the 'INCREASING KEY' or the 'DECREASING KEY' to set the demanded working temperature.

Finally, press again the ' SET' key for 5 seconds to confirm the setting parameters.

For example, the demanded working temperature is 37.0°C, and the former setting temperature is only 26.5°C. First, press the function key 'SET' on the control panel to enter the temperature setting program, till the setting temperature display SV is flashing. Then, press the 'SHIFTING KEY' to move the cursor to the second digital position before the base point, then press the 'INCREASING KEY' to change the display from '2' to '3'. Afterwards, again work with the 'SHIFTING KEY' and the 'INCREASING KEY' or the 'DECREASING KEY' to set the each digital position to make the display as 37.0°C as demanded working temperature. Finally, press again the 'SET' key for 5 seconds to confirm the setting parameters.

All the following functions are pre-set by default. They can only be modified after cancelling locking out of parameter modification = “LCK” function.

3. Upper limit temperature alarm setting

Before delivery, the upper limit temperature alarm setting of all units were set at 10°C. Therefore, normally it is not necessary to re-set it again. If it is necessary to re-set it again, then firstly press the function key 'SET' for 5 seconds, the controller enters the upper limit temperature setting program AU. Press the 'SHIFTING KEY' and the 'INCREASING KEY' or the 'DECREASING KEY' to set the demanded upper limit temperature alarm setting. Finally, press again the ' SET' key for 5 seconds to confirm the setting parameters.

4. Temperature display correction

Before delivery, all units were carefully tested ; normally it is not necessary to correct the temperature display. If the working environment is not good (such as the environment temperature is too high or too low), which makes the actual temperature display be different from the actual temperature, then it is necessary to correct it. The detailed steps are as following : firstly press the function key 'SET' for 5 seconds to enter the parameter setting cycle AL1. Then keep on pressing the 'SET' key, till the PV (Parameter Value) display 'SC' as temperature display correction. Then press the 'SHIFTING KEY' and the ' INCREASING KEY' or the 'DECREASING KEY' to set actual temperature display. Finally, press again the ' SET' key for 5 seconds to confirm the setting parameters.

5. After setting, all parameters will be stored for a long period of time. The unit will start heating up ; the heating indicator light will turn on. While closing to the setting temperature, the heating indicator light keeps on flashing, till the unit work at the setting temperature.

6. Open inner the door and put the things which must be cultured into the incubator, and the close them well. If the inner door stays open too long, it is a normal phenomenon that the incubator temperature will somewhat fluctuate.

7. After culture is finished, push the power switch to the '0' position. Please, don't open the incubator's door if the things are not to be taken out immediately.

8. If there is a high demand to temperature control accuracy and the temperature deviation, then the self-adjust PID function will be a help. When the temperature of the units is closed to the setting temperature, firstly press the function key, 'SET' for 5 seconds to enter the parameter setting cycle 'AL1' . Then keep on pressing the 'SET' key, till the PV display 'ATU' . The SV display '0000' , then press the 'INCREASING KEY' to make SV display as '0001'. Finally press again the 'SET' key for 5 seconds to confirm the setting parameters. Then the self adjust indicator light will turn on. The unit will enter the self-adjust PID control.

9. Heat power OUT- setting

Before product leave factory, the power was set the 50 %. In a general way, do not need set it again. If the environment temperature is too low or too high, which willlead to the temperature rises slowly. Than you can change the power to adjust. The range of power is 0 to 100% .Put "SET" key 8 seconds, then PV shows "HL" , put "SET" key again, PV shows "OUT" and SV shows "0050" .And then put the " < " , put " Δ " and " V " assorted operation. Lastly, put "SET" to sure it, the power-setting end.

10. Unlocking to modify the pre-settings :

Press the "SET" button for 5 seconds, the "AL1" menu will be displayed, press "SET" by scrolling the various functions until "LCK". Using the buttons "Δ" or "∇" displayed, correct the value to "0". Press "SET" to confirm. Locking out is now cancelled. The functions can be modified. To relock the commands and so avoid any accidental changes ; repeat the above procedure and select a value other than "0" for the "LCK" function.

7 Please pay attention that


1. The outer casing of the incubator must be effectively grounded so as to ensure safety in use.

2. The incubator should be placed in a room with well-ventilated conditions. Never put any inflammables or explosives around.


3. The incubator should never be too crowded with things to be cultured. There must be space in between.

4. The interior and exterior of the incubator should be always kept clean and tidy. If the incubator will not be used for long, it must be covered with a plastic dust cover and placed in a dry room.

8 Block-diagram of the electrical principle

9 Wiring diagram

10 Trouble shooting

Failure	Cause	Shooting method
1. No power supply	1. The plug has not been well inserted in the outlet, or the electric wire is broken.	1. Inserted the plug well in the outlet or connect the broken wire.
	2. Fuse is blown.	2. Replace the fuse.
2. The in-incubator temperature will no rise	1. The preset temperature is too low.	1. Adjust the preset temperature.
	2. The electric heater is worn out.	2. Replace the electric heater.
	3. The temperature controller is worn out.	3. Replace the temperature controller.
	4. The circulating fan is worn out.	4. Replace the fan.
3. A great deviation of the in-incubator temperature from the preset temperature	1. The temperature sensor is worn out.	1. Replace the temperature sensor.
	2. The trimmer potentiometer has not been well adjusted.	2. Adjust well the trimmer potentiometer.
4. Abnormal over temperature alarm	1. The preset temperature is too low.	1. Adjust the preset temperature.
	2. The temperature controller is worn out.	2. Replace the temperature controller.

Laboratory
Incubator 20 litres
Ref :
535 014

11 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact :

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