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
# 1 Safety

The equipment comes complete with safety features.

This manual indicates areas of possible risk.

## 2.1 Safety Icons

They identify risky situations and safety measures that should be adhered to. The icons make reference within the paragraph, marked as a lined grey.

### **Risk danger**

Risk danger

Beware, follow instructions as indicated.

### **Electrical risk**

Risk of electrical shock if accessing zones are shown with this icon.

Beware, follow instructions as indicated.

### **Risk of burns in high temperature areas.**

The temperature in the zone indicated with this icon can exceed 60°C. Use thermal protective gloves to complete described functions.

Beware, follow instructions as indicated

### **Important information**

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and maintain its optimum performance.

## 2.2 Risks that the operator can be subjected to

- Possibility of touching surfaces that are over 60°C.
- Possible escape of steam vapours.
- Electrical risk.

## 2.3 Qualified users

This equipment should only be used by personnel that are qualified users.

This equipment should only be used by personnel that have understood this manual or have previous adequate instruction of the use of steam sterilizers.

## 2 General information

2.1 Handle the equipment with care. Unpack and check that all items coincide with the delivery note. If you see any discrepancy or damage of any kind, notify your distributor as soon as possible.

2.2 Read this manual before using the equipment.

2.3 These instructions are an important part of the equipment and should be made available to all users.

2.4 If in any doubt or clarity on the use of this instrument, contact your local distributor or the service department of J.P. Selecta, s.a.u.

**2.5 ATTENTION!! EQUIPMENT NOT CLEAN AND DISINFECTED SHALL NOT BE REPAIRED.**

2.6 Any unauthorized modification, lack of maintenance or deliberate damage to the machine, transgresses the directive 89/655/CEE where the manufacturer will not be responsible for damages that could occur due to an undue care and attention.

2.7 Do not use the equipment for fluids that can produce vapours and that can form inflammable or explosive mixtures.

**DO NOT OPEN THE LID OF THE AUTOCLAVE IF THE MACHINE IS NOT CONNECTED TO THE MAINS SUPPLY AND THE MAINS SWITCH IS ILLUMINATED. IF NOT CONNECTED, DO NOT FORCE THE DOOR HANDLE.**

## 3 Technical specifications.

**Voltage supply 115-230V 50/60Hz according to the characteristics plate indications.**

Código Code	Capacidad L Capacity L	Medidas útiles CM Chamber dimensions CM		Medidas exteriores CM External dimensions CM			Rango °C Range °C	Potencia W Power W	Peso Kg Weight Kg
		Dimen.	Alto Height	Alto Height	Ancho Length	Fondo Width			
4002136	80L	40	60	107	58	72	115 a/to 134	3600	100
4002137	150L	50	70	118	80	95		7500	250

## 4 Contents list and accessories

### 4.1 Standard contents

The standard equipment consists of the following items:

5600-5603

34100

Descripción Description	Código Code	
Presoclave II	4002136	4002137
Cubrerresistencias Shelves	5600	5603
Racord salida vapor Steam outlet fitting	34100	
Manual instrucciones Instructions manual	80219	

1000496

### 4.2 Accessories

- Wired baskets made from AISI 304 stainless steel.

Código cestillo Basket code	Medidas CM Dimensions CM		Para código: To code:	Capacidad cestillos Baskets capacity
	Dimen.	Altura Height		
1000496	38	28	4002136	2
1000780	44	22	4002518	3

1001222

- Basket dividers made from AISI 304 stainless steel
  - 1001218 Divider 6 compartments for basket 1000496
  - 1001222 Divider 7 compartments for basket 1000780

## **5 Introduction**

The Presoclave II autoclaves are versatile and suitable for a wide range of applications, sanitation, industrial processes and quality control. These autoclaves can sterilize solids, with or without being wrapped, liquids, culture media.

### **5.1 Definitions**

#### **5.1.1 Sterilization**

Sterilization is the destruction or elimination of all forms of microorganisms, including spores present in inanimate objects.

#### **5.1.2 Disinfect**

Disinfection is the process of the destruction of microorganism agents. The temperature used is less than that of sterilization.

#### **5.1.3 Culture media**

A growth medium or culture medium is any substance in which microorganisms or cells can grow and multiply. Microorganisms can therefore be isolated and used in susceptibility testing. Generally produced as dry powder or granular form, but can also be produced as hydrates and preprepared. Should not be used unless they have been sterilized.

#### **5.1.4 Purge**

The purge cycle is when air is eliminated from the interior of the vessel with the objective of producing saturated steam vapour.

### **5.7 Saturated steam**

Water vapour at a temperature corresponding to the ebullition point of the liquid origin.

### **5.8 Atmospheric purge**

The term atmospheric purge is when air is expelled gravimetrically by vapour for a determined period.

## 6 Equipment description

The Presoclave II autoclave belongs to a series of saturated vapour autoclaves that are distinguished for the following features:


1. 6 point radial locking system.
2. Temperature and time electronic control circuit.
3. Safety blocking lock that only allows the opening of the lid with pressure 0 and temperature under 95°C.
4. Maximum pressure limiter that disconnects the heating element in cases of over-pressure.
5. Limiter thermostat to protect the heating element.

### 6.1 Equipment

- Safety system:
  - Safety overpressure valve.
  - Safety pressure switch.
  - Safety thermostat.
  - Lid closed detector.
  - Protective thermal lid cover.
- Chamber and lid made from stainless steel.
- External case made from stainless steel.
- Manual release vapour valve.
- Manual emptying valve.
- Waste water tank connection.

### 6.2 Features

- Sterilization temperature controllable from 115°C to 134°C.
- Temperature resolution 1°C.
- Sterilization time from 3 minutes to 99 minutes.
- Time resolution 1 minute.
- Suitable to work with solids and liquids atmospheric purge from 115°C to 134°C.

### 6.3 Control panel

1. Pressure gauge
2. Open lid indicator lamp
3. Open valve indicator lamp
4. End cycle indicator lamp
5. Overtemp alarm indicator lamp
6. Statement indicator lamp (see section 8.5 operation)
7. Temperature and errors screen
8. Time and pressure 0 screen
9. Temperature set key
10. Time set key
11. Increase value key
12. Decrease value key
13. Stop key
14. Starting up key

### 6.4 Elements position

15. Emptying manual valve control
16. Steam manual valve control
17. Open lid handle

18. Main switch

15. Emptying valve handle

16. Steam valve handle

19. Safety valve

20. Exit steam valve

21. Power supply cord

22. Characteristics plate indications

23. Thermostat reset button

24. Condenser outlet and chamber emptying

25. Open steam valve sensor

## 7 Installation

- 7.1 Locate the autoclave close to an appropriate power source.
- 7.2 The autoclave needs to be placed on a flat stable horizontal surface, leaving a minimum space of 10 cm around the instrument.
- 7.3 Adjust the front feet to level the autoclave and set the autoclave in position.

- 7.4 Fix a plastic or rubber hose to the pump output connection (24) with a receptacle to collect any condensate produced by the atmospheric purge.

**¡WARNING! DON'T TRY TO OPEN THE LID UNLESS READING THE «OPERATION» SECTION**

## 8 Operation

### 8.1 Opening the Lid

- The autoclave must be connected to mains supply and the main switch turned on.
- Pressure gauge: 0 bar
- Temperature in the internal chamber must be under 95°C.
- The steam manual valve (16) must be opened.

**ONLY WITH THIS CONDITIONS LID IS READY TO OPEN**



### 8.2 Start

- Start on using the main switch (18)
- Open the lid as we can see on 8.1 section, using the handle (17)
- Place the shelf at the bottom of chamber.
- Be sure the valve (15) is closed.
- Fill of decalcified water until shelf level.
- Put the material to sterilize and close the lid.
- Close the manual emptying valve (16)
- Start the cycle as shows the section 8.3

**BEFORE STARTING THE CYCLE, CHECK ALWAYS THERE IS ENOUGH WATER IN THE CHAMBER; IF NOT, THE HEATING ELEMENT COULD BE DAMAGED.**

**8.3 Temperature Cycle Selection**

- To select the working temperature (from 115 to 134°C), press «°C» key (9) and, simultaneously, the increase value key (11) or the decrease value key (12). Temperature is shown on display.

**8.4 Time Cycle Selection**

- To select the cycle time (from 3 min to 99 min), press the set time key (10) and, simultaneously, the increase value key (11) or the decrease value key (12). Time is shown on display.

**8.5 Start Working Cycle**

- Once time finished previous operations on sections 8.2, 8.3 and 8.4, press the start on key (14) and autoclave will work with the following indications:

Lamp (6): Lighting  
 Display (7): Real temperature is shown.  
 Display (8): 6 blinking horizontal lines.

**8.6 Working Cycle Steps**

- 1 Purge: The autoclave expulses the chamber internal air to get saturated steam. The Presoclave II autoclaves makes a process called atmospheric purge. The temperature increases up to reach 103°C. The purge electrovalve holds opened during 3 minutes once 103°C have been reached.

Lamp (6): Lighting  
 Display (7): Show temperature below 103°C.  
 Display (8): 6 blinking horizontal lines.

- 2 Temperature increasing until the set temperature: The temperature continues increasing from 103°C to the set temperature, with all the valves closed.

Lamp (6): Lighting  
 Display (7): Real temperature is shown between 103°C and set temperature.  
 Display (8): 6 blinking horizontal lines.

- 3 Sterilization: 20 seconds after reaching the set temperature the autoclave starts the sterilization during all the time that we have set.

Lamp (6): Permanent lighted on.  
 Display (7): Real temperature is shown, the same set temperature.  
 Display (8): Remain minutes to finish the cycle.

**DONOT OPENTHE VALVES 15 & 16 UNTIL THE CYCLE IS COMPLETELY FINISHED TO AVOID POSSIBLE BURNS WITH HEAT WATER OR STEAM.**

- End cycle: When sterilization is finished, cold process starts. The lamp (6) is blinking and the lamp (4) lights on. The display (7) shows the internal real temperature.  
In case of solids sterilization fast steam output can be made using the manual valve (16), the steam output is made with output (20)


## ¡WARNING!

- **THE DIRECT CONTACT WITH PRESSURE STEAM CAN PRODUCE SERIOUS BURNS.**
- **IF YOU USE A HOSE TO MAKE THE STEAM OUTPUT, DO NOT HOLD IT WITH THE HANDS TO AVOID BURNS.**

**YOU CAN NOT OPEN THE LID UNTIL THE PRESSURE IS NOT 0 BAR AND THE INTERNAL TEMPERATURE UNDER 95°C. SEE SECTION 8.1**

When the lid is ready to be opened the display (7) shows the word «End» and the lamp (6) is off.

If liquids or culture media are sterilized must leave cooler until the pressure be 0 bar. During the cooling, is the same case of manual depressurize, the display shows the internal temperature.

When the pressure is 0, the temperature is under 95°C, the display (7) shows the word «End». Therefore, the valve (20) can be opened by the handle (16) and open the lid with his handle (17).

## 9 Error messages

**Error 1** Pressure sensor failure or overpressure.

**Error 2** PT100 probe failure or overtemp.

**Error 3** External alarm. Pressure switch or safety thermostat triggered.

**To reset the safety thermostat press its button number 23**

The thermostat is adjusted in the factory it must not be operated except to reset it.

**Error 4** Manual valve (16) opened during a cycle.

**Error 5** Opened door during a cycle.

**Error 6** Pressure over +0.05bar with the autoclave stopped.

**Error 7** E2PROM failure.

**Error 8** Out of limits temperature during the sterilization.

**Error 9** Mains supply failure during cycle.

## 10 Spare parts

Presoclave II 80 litres. Part No. 4002136

07292	RS-232 communication wire
07315	26 pins communication wire
13006	Contactora
15411	Fuse 16 A
16055	Solenoid electromagnet
16057	Electrovalve NA 3 mm 24 V
16212	Pressure switch
20002	Lid microruptor
20067	Valve microruptor
20096	Bipolar switch
29358	Frontal circuit
29362	Power circuit
36027	Solid relay 25A
39001	Heating element 3680W/230V
43031	Thermostat
43071	Pt100 probe

Presoclave II 150 litres. Part No. 4002137

## 11 Recommendation to achieve a perfect sterilization

The material to be sterilized needs to be clean, free from any encrusted or residual material, we recommend that they are washed in a good detergent and distilled water, then rinsed with distilled water.

Do not overload the baskets, drums or trays, always leave a space of 1 to 2 cm between items to allow a free passage of vapour and to ease drying.

### 11.1 Sterilization of instruments in bulk

Place the instruments on a tray with water repellent paper leaving a space between each instrument.

Try to keep similar type material instruments together on a tray.

### 11.2 Tubing

The tubing need to be clean and drained. Place the tubing so that the openings on both ends are open. Make sure that the tubes are not touching the autoclave wall, blocked or twisted and that nothing can block the inside.

### 11.3 Vessels

Never put hermetically sealed vessels in the autoclave.

Place the contents open end down to prevent any residual water being left inside.

### 11.4 Liquids


Place the receptacle with the liquid inside on a drip tray, making sure that the receptacle can withstand the temperature.

Do not fill the receptacle more than 2/3rds of its maximum capacity, DO NOT CLOSE THE LID IF IT HAS A HERMETICAL SEALED, only use a lid if the pressure can be released. (eg. Serotap or Bacti-lid)

Try to avoid using narrow neck vessels.


At the end of the cycle, let the autoclave cool down slowly (until the pressure is 0 kg/cm<sup>2</sup>)

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